



Product Datasheet

PHOX2B antibody (orb395436)

Description	Mouse monoclonal antibody to PHOX2B.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	paired-like homeobox 2b
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	35 kDa
Uniprot ID	Q99453
Dilution Range	WB: 1:1000-1:5000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis
of Neuro-2a cells ...



Immunohistochemical
staining of paraffin...