


Product Datasheet

ZC3HC1 antibody (orb39544)

Description	ZC3HC1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Nuclear-interacting partner of ALK protein (1-260AA)
Target	ZC3HC1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86WB0
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
HeLa whole
cel...