



Product Datasheet

PI3K p55 (gamma) antibody (orb395442)

Description	Rabbit polyclonal antibody to PI3K p55 (gamma).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	phosphoinositide-3-kinase, regulatory subunit 3 (gamma)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	55 kDa
Uniprot ID	Q92569
Dilution Range	WB: 1:500 - 1:1000; IHC: 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tis...



Immunohistochemical staining of paraffin...