



Product Datasheet

PI3K p85 (alpha) antibody (orb395443)

Description	Mouse monoclonal antibody to PI3K p85 (alpha).
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	phosphoinositide-3-kinase, regulatory subunit 1(alpha)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	85 kDa
Uniprot ID	P27986
Dilution Range	WB:1:1000-1:10000, IP:1:200-1:2000, IHC:1:20-1:200, IF/ICC:1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Neuro-2a cells ...



Immunohistochemical staining of paraffin...