



Product Datasheet

PTPRCAP antibody (orb395488)

Description	Rabbit polyclonal antibody to PTPRCAP.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	protein tyrosine phosphatase, receptor type, C-associated protein
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	30 kDa
Uniprot ID	Q14761
Dilution Range	WB: 1:200-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells us...



Immunohistochemical staining of paraffin...