

Product Datasheet

RAB5A-Specific antibody (orb395499)

Description Rabbit polyclonal antibody to RAB5A-Specific.

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen RAB5A, member RAS oncogene family

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 24 kDa


Uniprot ID **P20339**

Dilution Range WB: 1:500-1:3000; IHC: 1:50-1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of human brain tis...



Immunohistochemical staining of paraffin...