


Product Datasheet

S1 tag antibody (orb395539)

Description	Mouse monoclonal antibody to S1 tag.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	peptide
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	68 kDa
Dilution Range	WB:1:5000-1:50000
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant pro...