



Product Datasheet

SERPINA10 antibody (orb395569)

Description	Mouse monoclonal antibody to SERPINA10.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	serpin peptidase inhibitor, clade A(alpha-1 antiproteinase, antitrypsin), member 10
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	51 kDa
Uniprot ID	Q9UK55
Dilution Range	WB: 1:1500-1:6000; IF: 1:50-1:500; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using SERPINA10 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using SERPINA10 antibody.