


Product Datasheet

SLA/LP antibody (orb395579)

Description	Rabbit polyclonal antibody to SLA/LP.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Sep(O-phosphoserine) tRNA:Sec(selenocysteine) tRNA synthase
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	56 kDa
Uniprot ID	Q9HD40
Dilution Range	WB: 1:500-1:2000; IP: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using...