

Product Datasheet

SMYD1-Specific antibody (orb395592)

Description Rabbit polyclonal antibody to SMYD1-Specific.

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IHC, IP, WB

Immunogen SET and MYND domain containing 1

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 57 kDa


Uniprot ID **Q8NB12**

Dilution Range WB: 1:500-1:1000; IP: 1:500-1:1000; IHC: 1:20-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of mouse heart tis...



Immunohistochemical staining of paraffin...