



Product Datasheet

SSR4 antibody (orb395613)

Description	Rabbit polyclonal antibody to SSR4.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	signal sequence receptor, delta(translocon-associated protein delta)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	19 kDa
Uniprot ID	P51571
Dilution Range	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using...



Immunohistochemical staining of paraffin...