



Product Datasheet

ST6GALNAC6 antibody (orb395614)

Description	Mouse monoclonal antibody to ST6GALNAC6.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	ST6(alpha-N-acetyl-neuraminy-2,3-beta-galactosyl-1,3)-N-acetylgalactosaminide alpha-2,6-sialyltransferase 6
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	40 kDa
Uniprot ID	Q969X2
Dilution Range	WB: 1:500-1:1000; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO 320 cells ...



Immunohistochemical staining of paraffin...