

## Product Datasheet

### Stanniocalcin 2 antibody (orb395617)

**Description** Mouse monoclonal antibody to Stanniocalcin 2.

**Species/Host** Mouse

**Reactivity** Human, Mouse

**Conjugation** Unconjugated

**Tested Applications** ELISA, IF, IHC, WB

**Immunogen** stanniocalcin 2

**Form/Appearance** liquid

**Storage** PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

**Note** For research use only

**Isotype** IgG2a

**Clonality** Monoclonal


**Purity** ≥95% as determined by SDS-PAGE

**MW** 33 kDa


**Uniprot ID** **O76061**

**Dilution Range** WB:1:500-1:5000, IHC:1:20-1:200, IF/ICC:1:10-1:100

**Expiration Date** 12 months from date of receipt.



Western blot analysis of HEK-293 cells u...



Immunohistochemical staining of paraffin...