



Product Datasheet

STT3B antibody (orb395625)

Description	Rabbit polyclonal antibody to STT3B.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	STT3, subunit of the oligosaccharyltransferase complex, homolog B(<i>S. cerevisiae</i>)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	94 kDa
Uniprot ID	Q8TCJ2
Dilution Range	WB: 1:1000-1:4000; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain tis...



Immunohistochemical staining of paraffin...