



Product Datasheet

SULT1A4 antibody (orb395629)

Description	Rabbit polyclonal antibody to SULT1A4.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	sulfotransferase family, cytosolic, 1A, phenol-preferring, member 4
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	32 kDa
Uniprot ID	PODMNO
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of human colon tis...



Immunohistochemical staining of paraffin...