



Product Datasheet

Syntaxin 6 antibody (orb395655)

Description	Mouse monoclonal antibody to Syntaxin 6.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	syntaxin 6
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	29kd
Uniprot ID	043752
Dilution Range	WB: 1:2000-1:15000; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells us...



Immunohistochemical staining of paraffin...