

Product Datasheet

TAF1-Specific antibody (orb395660)

Description Rabbit polyclonal antibody to TAF1-Specific.

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen TAF1 RNA polymerase II, TATA box binding protein(TBP)-associated factor, 250kDa

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 50 kDa


Uniprot ID **P21675**

Dilution Range WB: 1:500-1:3000; IHC: 1:20-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of mouse brain tis...



Immunohistochemical staining of paraffin...