



Product Datasheet

TBC1D5 antibody (orb395662)

Description	Rabbit polyclonal antibody to TBC1D5.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	TBC1 domain family, member 5
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	89 kDa
Uniprot ID	Q92609
Dilution Range	WB: 1:500-1:2400; IP: 1:200-1:1000 ; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 cells using TBC1D5 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using TBC1D5 antibody.