



Product Datasheet

TLE3 antibody (orb395681)

Description	Mouse monoclonal antibody to TLE3.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	transducin-like enhancer of split 3(E(sp1) homolog, Drosophila)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	90 kDa
Uniprot ID	Q04726
Dilution Range	WB: 1:500-1:2000; IP: 1:500-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK-293 cells u...



Immunohistochemical staining of paraffin...