



Product Datasheet

TMEM195 antibody (orb395684)

Description	Rabbit polyclonal antibody to TMEM195.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	transmembrane protein 195
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	52 kDa
Uniprot ID	Q6ZNB7
Dilution Range	WB: 1:500-1:1000; IHC: 1:100-1:400
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 cells using TMEM195 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using TMEM195 antibody.