



Product Datasheet

TMF1-Specific antibody (orb395690)

Description	Rabbit polyclonal antibody to TMF1-Specific.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	TATA element modulatory factor 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	150 kDa
Uniprot ID	P82094
Dilution Range	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using TMF1-Specific antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using TMF1-Specific antibody.