



Product Datasheet

TOMM7 antibody (orb395696)

Description	Rabbit polyclonal antibody to TOMM7.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	translocase of outer mitochondrial membrane 7 homolog(yeast)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	6 kDa
Uniprot ID	Q9P0U1
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis
of Raji cells using...



Immunohistochemical
staining of paraffin...