


Product Datasheet

TPI1 antibody (orb395698)

Description	Mouse monoclonal antibody to TPI1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	triosephosphate isomerase 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	27 kDa
Uniprot ID	P60174
Dilution Range	WB: 1:1000-1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using TPI1 antibody (orb395698).