


Product Datasheet

transcription termination factor-Specific antibody (orb395703)

Description	Mouse monoclonal antibody to transcription termination factor-Specific.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	transcription termination factor, RNA polymerase I
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	43 kDa, observed 43 kDa, 66 kDa
Uniprot ID	Q15361
Dilution Range	WB:1:1000-1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using the transcription termination factor-Specific antibody.