



Product Datasheet

Transgelin-2-specific antibody (orb395706)

Description	Rabbit polyclonal antibody to Transgelin-2-specific.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, IP, WB
Immunogen	transgelin 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	22 kDa
Uniprot ID	P37802
Dilution Range	WB:1:200-1:2000, IP:1:200-1:2000, IHC:1:20-1:200, IF/ICC:1:10-1:100, FACS:N/A
Expiration Date	12 months from date of receipt.



Western blot analysis of rat stomach tis...



Immunohistochemical staining of paraffin...