



Product Datasheet

TRAPB antibody (orb395711)

Description	Rabbit polyclonal antibody to TRAPB.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	signal sequence receptor, beta(translocon-associated protein beta)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	34 kDa
Uniprot ID	P43308
Dilution Range	WB:1:500-1:5000, IHC:1:20-1:200, IF/ICC:1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using TRAPB antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using TRAPB antibody.