


Product Datasheet

TRBV5-4 antibody (orb395714)

Description	Rabbit polyclonal antibody to TRBV5-4.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	T cell receptor beta variable 5-4
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	39 kDa, observed 45 kDa
Uniprot ID	A0A0C4DH59
Dilution Range	WB:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of U-937 cells using...