



Product Datasheet

TRPC4AP antibody (orb395717)

Description	Mouse monoclonal antibody to TRPC4AP.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	transient receptor potential cation channel, subfamily C, member 4 associated protein
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	91 kDa
Uniprot ID	Q8TEL6
Dilution Range	WB:1:500-1:5000, IHC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis
of human heart tis...



Immunohistochemical
staining of paraffin...