



Product Datasheet

TSG-6 antibody (orb395718)

Description	Rabbit polyclonal antibody to TSG-6.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	tumor necrosis factor, alpha-induced protein 6
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	41 kDa
Uniprot ID	P98066
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH/3T3 cells u...



Immunohistochemical staining of paraffin...