



Product Datasheet

TSNAXIP1 antibody (orb395720)

Description	Rabbit polyclonal antibody to TSNAXIP1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	translin-associated factor X interacting protein 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	33 kDa
Uniprot ID	Q2TAA8
Dilution Range	WB:1:200-1:2000, IHC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse testis ti...



Immunohistochemical staining of paraffin...