



Product Datasheet

TTC37 antibody (orb395724)


Description	Rabbit polyclonal antibody to TTC37.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	tetratricopeptide repeat domain 37
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	35 kDa
Uniprot ID	Q6PGP7
Dilution Range	WB:1:200-1:1000, IHC:1:20-1:200, IF/ICC:1:10-1:100, IP:1:200-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using TTC37 antibody.



Immunohistochemical staining of human tissue using TTC37 antibody.



Immunoprecipitation analysis of mouse tissue using TTC37 antibody.