



Product Datasheet

TYW1B antibody (orb395735)

Description	Rabbit polyclonal antibody to TYW1B.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	tRNA-yW synthesizing protein 1 homolog B(non-protein coding)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	75 kDa
Uniprot ID	Q6NUM6
Dilution Range	WB: 1:200-1:2000, IHC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver tis...



Immunohistochemical staining of paraffin...