




## Product Datasheet

### TP73 antibody (orb39576)

|                            |                                                                                                                                                                                                                                                                                                                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Peptide affinity purified rabbit polyclonal antibody.                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Tested Applications</b> | FC, IF, IHC-P, WB                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Target</b>              | TP73                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.                                                                                                                                                                                                                                          |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>MW</b>                  | 69623                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Uniprot ID</b>          | <a href="#">O15350</a>                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>NCBI</b>                | <a href="#">NP_001191121.1</a> , <a href="#">NP_001191119.1</a> ,<br><a href="#">NP_001119713.1</a> , <a href="#">NP_001191114.1</a> ,<br><a href="#">NP_001191118.1</a> , <a href="#">NP_001119712.1</a> ,<br><a href="#">NP_001191113.1</a> , <a href="#">NP_001191115.1</a> ,<br><a href="#">NP_005418.1</a> , <a href="#">NP_001119714.1</a> ,<br><a href="#">NP_001191120.1</a> , <a href="#">NP_001191117.1</a> |
| <b>Dilution Range</b>      | WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50, IF/ICC: 1:10-50                                                                                                                                                                                                                                                                                                                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                                                                                                                                                                       |



Immunofluorescence analysis of 293 cell ...



Immunohistochemical staining of paraffin...