


Product Datasheet

XLF antibody (orb395773)

Description	Rabbit polyclonal antibody to XLF.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	nonhomologous end-joining factor 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	33 kDa
Uniprot ID	Q9H9Q4
Dilution Range	WB:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of human small int...