



Product Datasheet

ZGLP1 antibody (orb395775)

Description	Rabbit polyclonal antibody to ZGLP1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	GATA like protein-1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	30 kDa
Uniprot ID	P0C6A0
Dilution Range	WB: 1:200-1:2000, IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse ovary tis...



Immunohistochemical staining of paraffin...