


Product Datasheet

ZIP kinase antibody (orb395778)

Description	Rabbit polyclonal antibody to ZIP kinase.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	death-associated protein kinase 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	68 kDa
Uniprot ID	O43293
Dilution Range	IF/ICC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK-293 cells u...