



Product Datasheet

ZIP4 antibody (orb395779)

Description	Rabbit polyclonal antibody to ZIP4.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	solute carrier family 39(zinc transporter), member 4
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	50kd,56kd
Uniprot ID	Q6P5W5
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK-293 cells u...



Immunohistochemical staining of paraffin...