



Product Datasheet

ATRX antibody (orb395790)


Description	ATRX antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human ATRX.Å
Target	ATRX
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	282 kDa
Uniprot ID	P70486, P46100, Q61687
NCBI	000480, 000489
Entrez	546
Dilution Range	IHC-P: 1:50-300
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



IHC-P staining of mouse liver tissue
usi...



Immunohistochemical staining of human
lu...