

Product Datasheet

Wnt-5a antibody (orb395793)

Description

Wnt-5a antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Wnt-5a. Please contact us for the exact immunogen sequence.

Target

Wnt-5a

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml & - 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

42 kDa

Uniprot ID

[P41221](#), [P22725](#), [Q9QXQ7](#)

NCBI

[001243034](#), [001256105](#)

Entrez



[7474](#)

Dilution Range


IHC-P: 1:50-300, WB: 1:100-800

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of rat thyr...



WB analysis of mouse lung (lane 1), rat ...



IHC-P image of rat thyroid tissue using

...