


Product Datasheet

KPNA2 antibody (orb395799)

Description	KPNA2 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human KPNA2 . Please contact us for the exact immunogen sequence.
Target	KPNA2
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	58 kDa
Uniprot ID	P52293, P52292, Q6P6T9
Entrez	3838
Dilution Range	WB: 1:200-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat lung (lane ...)