


Product Datasheet

LYVE1 antibody (orb395802)

Description	LYVE1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human LYVE1. Please contact us for the exact immunogen sequence.
Target	LYVE1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	35 kDa
Uniprot ID	Q8BHC0 , Q9Y5Y7 , D3ZD19
NCBI	006682 , 006691
Entrez	10894
Expiration Date	12 months from date of receipt.



WB analysis of rat heart (lane 1), human...