



Product Datasheet

PTEN antibody (orb395804)

Description	PTEN antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PTEN. Please contact us for the exact immunogen sequence.
Target	PTEN
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	47 kDa
Uniprot ID	054857 , P60484
NCBI	000314 , 000305
Entrez	5728
Dilution Range	IHC-P: 1:50-300, WB: 1:100-1000
Expiration Date	12 months from date of receipt.



WB analysis of rat stomach (lane 1), mou...



IHC-P image of rat stomach tissue using ...