




Product Datasheet

AIFM1 antibody (orb39603)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	AIFM1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	66901
Uniprot ID	O95831
NCBI	NP_001124318.1 , NP_665811.1 , NP_001124319.1 , NP_004199.1 , NP_665812.1
Dilution Range	IF/ICC: 1:200, WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U251 cell...



Immunofluorescence analysis of U251 cell...