

Product Datasheet

ANKRD36B antibody (orb39612)

Description	ANKRD36B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Ankyrin repeat domain-containing protein 36B protein (1-230AA)
Target	ANKRD36B
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N2N9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.