

Product Datasheet

IL-6 receptor antibody (orb396139)

Description

Mouse monoclonal antibody to IL-6 receptor.

Species/Host

Human

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

The original antibody was generated by CDR grafting of mouse PM-1 antibody onto human REI and NEW framework regions to create rhPM-1.

Target

IL-6 receptor

Preservatives

PBS with 0.02% Proclin 300.

Concentration

1 mg/ml

Storage

Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.

Note

For research use only

Application notes

This antibody is derived from the parental mouse clone PM-1 (PM1). Tocilizumab (TCZ) binds and neutralizes IL-6R, resulting in the inhibition of various IL-6-mediated biological activities, including inflammation-related, immunomodulatory and tissue/matrix remodelling effects. This antibody is extensively explored for its therapeutic potential in treatment of rheumatoid arthritis (PMID: 25099958; 9836380). Tocilizumab can be used in ELISA assays and is able to inhibit the growth of IL-6-dependent tumor cell lines. It inhibits sIL-6R and mIL-6R-mediated signaling. A study reported that a combination of tocilizumab and methotrexate (MTX) results in greater inhibition of joint damage and improvement in physical function than does MTX alone (PMID: 21360490). Treatment with tocilizumab results in significant improvement, maintained over time, of polyarticular-course juvenile idiopathic arthritis (pcJIA) signs and symptoms and has a safety profile consistent with that for adults with rheumatoid arthritis (PMID: 24834925).

Isotype

IgG1

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558