

Product Datasheet

DDDDK-tag antibody (orb396143)

Description	Mouse monoclonal antibody to DDDDK-tag.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Interleukin 2 genetically fused to the peptided DYKDDDDK, corresponding to the DDDDK epitope tag.
Target	DDDDK-tag
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody recognises the widely used DDDDK-tag commonly used as an epitope tag in fusion proteins. The tag can be recognised at any position within a fusion protein, i.e. N-terminally, Met-N-terminally, C-terminally as well as internally. Binding of this antibody is not sensitive to calcium. This antibody recognises both native and denatured protein.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.