

Product Datasheet

ERp57 antibody [Out of stock] (orb396317)

Description

Mouse monoclonal antibody against

Species/Host	Mouse
Reactivity	Human
Conjugation	ATTO-390
Tested Applications	ICC, WB
Immunogen	Human recombinant ERp57 (Grp58) protein corresponding to amino acids 25 - 505 (his-tagged)
Target	ERp57
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-266 was sufficient for detection of ERp57 in 15 µg of HeLa by ECL immunoblot analysis using Goat Anti-Mouse IgG:HRP as the secondary Antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	57 kDa
Uniprot ID	P30101
NCBI	NP_005304.3
Entrez	2923
Dilution Range	WB (1:1000); ICC/IF (1:100)
Expiration Date	12 months from date of receipt.