

Product Datasheet

CD74 antibody (orb396345)

Description

Mouse monoclonal antibody against CD74

Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ICC, WB
Immunogen	Human CD74 invariant chain synthetic peptide.
Target	CD74
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-267 was sufficient for detection of CD74 in 15 µg of B-Cell Lymphoma cell line Raji by ECL immunoblot analysis using Goat Anti-Mouse IgG:HRP as the secondary Antibody.
Isotype	IgG2a
Clonality	Monoclonal
MW	33 kDa
Uniprot ID	P04233
NCBI	NP_001020329.1
Entrez	972
Dilution Range	WB (1:1000); ICC/IF (1:100)
Expiration Date	12 months from date of receipt.