



Product Datasheet

Citrulline antibody (orb396417)


Description	Mouse monoclonal antibody against Citrulline conjugated to Biotin.
Species/Host	Mouse
Reactivity	All
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Synthetic L-Citrulline conjugated to Keyhole Limpet Hemocyanin (KLH).
Target	Citrulline
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-501 was sufficient for detection of Citrulline in 2 µg of Citrulline conjugated to BSA by ECL immunoblot analysis using Goat Anti-Mouse IgG:HRP as the secondary Antibody.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB (1:1000); ICC/IF (1:50); ELISA (1:1000); IHC: (1:25) optimal dilutions for assays should be determined by the user.
Expiration Date	12 months from date of receipt.



Immunohistochemistry analysis of mouse C...



Western Blot analysis of Lane 1: Molecul...



Immunofluorescence analysis of Embryonic...