



Product Datasheet

TGN46 antibody (orb39649)

Description	Rabbit polyclonal antibody to TGN46
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	TGOLN2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	45880
Uniprot ID	O43493
NCBI	NP_001193769.1 , NP_001193773.1 , NP_006455.2 , NP_001193770.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
MDA-MB231
ce...